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

RAPID SARS-COV-2 ANTIGEN TEST CARD  
FOR COVID-19 ANTIGEN SELF TESTING USE

Convenient. Easy to use. Results in 20 minutes

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
## Test Procedure


### For Nasal Swab Specimen


1. Open the Sample extraction buffer
2. Add all the buffer to the Extraction tube
3. Place the swab with specimen into the extraction tube.
4. Leave the swab in the extraction buffer for 1 minute.
5. Cover the extraction tube
6. Add 2-3 drops into the sample hole (S).
7. Read the result at 15-20 minutes.

### For Sputum Specimen

1. Open the cover of the dalusur tube.
2. Scrow on the salivary fummal.
3. Mzknz Kimai speed in the throat to clear alliva from to throat
4. Collect the Halita to 2 ml.
5. Scriw off the alitary funnel.
6. Cover and turs up tide down and mix wall
7. Screw off the cover.suck a tube of liquid with a dropper
8. Drop 2-3 drops into the sample hole and start to count for 10-15 minutes.

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**CE** 1434

## DECLARATION OF CONFORMITY

**According Directive 98/79/EC on in vitro diagnostic medical devices.**

**Manufacturer:** BioDetect (Xiamen) Biotechnology Co., Ltd.

**Address:** Building B9, Xiamen Biomedical Industrial Park, 2066 Wengjiao West Road, 361026, Xiamen, Fujian, China.

**European Representative:** Lotus NL B.V.

**Contact person:** Peter E-mail: peter@lotusnl.com

**Address:** Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

**In Vitro Diagnostic Directive:**

**Rapid SARS-CoV-2 Antigen Test Card**

**Category:** Self-Testing.

**Model:** 1 test/kit, 5 tests/kit, 20 tests/kit

**Conformity assessment route:** Declaration of Conformity IVDD Annex III.

**General applicable directive:**

Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices.

Notified Body: Polish Centre for Testing and Certification

Identification number: CE1434

EC Certificate No.: 1434-IVDD-521/2021

Start of CE Marking: 2021-12-22

We, the manufacturer, here declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In Vitro Diagnostic Medical Devices.

We agree to develop, implement and maintain a documented post-production monitoring process.

Represented by

Signed on:   
Place: Xiamen, China

BioDetect (Xiamen) Biotechnology Co., Ltd.

**Name of authorized signatory:** Li Boan

**Position held in the company:** General Manager

# EC Certificate

## CERTIFICATE

**EC Certificate No. 1434-IVDD-521/2021**

EC Design-examination  
Directive 98/79/EC concerning  
*in vitro* diagnostic medical devices

Polish Centre for Testing and Certification certifies  
that manufactured by:

**BioDetect (Xiamen) Biotechnology Co.,Ltd.**  
**Building B9 and B5, Xiamen Biomedical Industrial Park,**  
**2066 Wengjiao West Road, Xiamen, Fujian, P.R.China**

*in vitro* diagnostic medical devices  
self-testing

### RAPID SARS-COV-2 ANTIGEN TEST CARD

in terms of design documentation, comply with requirements  
of Annex III (Section 6) to Directive 98/79/EC (as amended)  
implemented into Polish law,  
as evidenced by the audit conducted by the PCBC  
Validity of the Certificate: from 21.12.2021 to 27.05.2024  
The date of issue of the Certificate: 21.12.2021  
The date of the first issue of the Certificate: 21.12.2021

Issued under the Contract No. MD-112/2021  
Application No: 169/2021  
Certificate bears the qualified signature.  
Warsaw, 20.12.2021  
Module A1  
FBM-30-E\_10

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Director  
Medical Device Certification  
Department

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
# Malaysia Registration

Official Portal  
 Medical Device Authority  
 Ministry of Health Malaysia

Medical Device Authority  
 Ministry of Health Malaysia

Our Hotline +603 - 8230 0300

Home Introduction Division FAQ Contact Us Documents Counter Appointment search...

Home / Announcement / SELF-TEST COVID-19 TEST KIT FOR CONDITIONAL APPROVAL (APPROVED)

## SELF-TEST COVID-19 TEST KIT FOR CONDITIONAL APPROVAL (APPROVED)

### SELF-TEST COVID-19 TEST KIT FOR CONDITIONAL APPROVAL (APPROVED)

List of RTK ANTIGEN (Self-Test) that is approved for Conditional Approval based on the decision on the consensus of the Covid-19 Test Kit Expert committee is as follows.

All test submissions are scored according to: -

- the manufacturer reported clinical and analytical performance evidence.
- the evaluation test results from testing facilities are according to the committee evaluation criteria set by Clinical expert panels.
- Supporting Documents for COVID-19 IVD Test Kits Conditional Approval.

The use of COVID-19 self- test kit shall be limited for screening purpose only and all test results need further confirmation using RT-PCR.

Below is the list of all tests that have been selected to date and the status is Conditional Approval (**please note: list is updated on a routine basis**).

Notes:  
 \*Sample type is based on testing facility evaluation report.

| NO | COMPANY NAME        | PRODUCT NAME  | MANUFACTURER  | IDENTIFIER | LOT NUMBER |
|----|---------------------|---|---|------------|------------|
| 79 | Dyna Source Sdn Bhd | BioDetect Rapid SARS-CoV-2 Antigen Test Card (Self-test)<br>Sample type: Nasal swab | BioDetect (Xiamen) Biotechnology Co., Ltd.<br>P.R.China | GMT102     | T220050901 |


## Production workshop

Our 100,000-class clean and dust-free workshop strictly implements ISO13485:2016 and GMP requirements. Fully automated production line, high-efficiency production technology, and strict quality control. The current daily production capacity can reach 2 million to 3 million.

## RAPID SARS-CoV-2 Antigen Test Card

| Product Picture   | Unit         | Specification | Dimensions        | Gross weight |
|---|--------------|---------------|-------------------|--------------|
|  | 1<br>Package | 1 Test        | 190mm*50mm*24mm   | 26.6g        |
|  | 1 Carton     | 250 Tests     | 648mm*394mm*274mm | 7.3kg        |

| Product Picture   | Unit         | Specification | Dimensions        | Gross weight |
|---|--------------|---------------|-------------------|--------------|
|   | 1<br>Package | 5 Tests       | 195mm*105mm*30mm  | 85.4g        |
|  | 1 Carton     | 500 Tests     | 660mm*410mm*230mm | 9.8kg        |